

Home Assignment Week-8 Day-1

Test case Generation with RAG

Target User Story / Query Specification

Select a sample teaching story or type your own custom user story.

Quick Presets:

WhatsApp Diagnostic Reports (HC-125)

Patient UHID Registration (HC-089)

Lab Report Verification & Alerting (HC-305)

Jira URL

Email

API Token

Project Key

Test connection

Fetch stories

Fetch & Ingest Project

https://ninahire.atlassian

ninahire@outlook.com

.....

KAN

Filter stories

Stories

Ingest

Auto Use

Yes

Select a story to preview its details here.

Module: Patient Communication & Diagnostics

User Story ID: HC-125

User Story Title: Share Diagnostic Reports with Patients via WhatsApp

User Story Description:

As a diagnostic technician, I want to automatically send patients their lab test reports securely through WhatsApp once results are validated, so that patients can access their diagnostic data conveniently without logging into the hospital portal.

Generated Test Cases with RAG Pipeline

Mode: With RAG

Model: openai/gpt-oss-120b

50 MongoDB Matches

3,714 / 8000 Tok (46% Budget Safe)

Export CSV

Generated Test Cases (2)

Retrieved MongoDB Context (8)

Linked User Stories (10)

ID	Title	Module	Priority	Type	Preconditions	Test Steps	Expected Results
TC_6355	Successful delivery of encrypted diagnostic PDF via WhatsApp after validation	Patient Communication & Diagnostics	P1	Integration	Technician - • Patient record contains a valid, consent-checked WhatsApp number (OTP-verified). • Lab report status = "Validated". • WhatsApp API service is reachable and health-checked. • HIPAA-compliant encryption library is configured.	4. Click 'Send Report via WhatsApp'. 5. System generates an encrypted PDF (AES-256) and creates a secure download URL with a 24-hour TTL. 6. System invokes WhatsApp Business API to send a templated message containing the URL to the patient's verified mobile number. 7. Verify that the message status returned by the API is 'sent'. 8. From a separate device, click the URL within 23 hours and confirm the PDF downloads, prompts for password, and decrypts correctly. 1. Open Laboratory → Post Analytical Phase → Report Result for patient UHID 789012.	Secure URL contains a UUID token and expires exactly after 24 hours. • WhatsApp API response = success (message_id returned). • Audit log entry created with fields: UHID, PRN, Report ID, delivery method=WhatsApp, timestamp, status=Sent. • Patient can download and open the PDF before expiration; after 24 hours the URL returns HTTP 410 Gone. • Consent check

Test case Generation without RAG

Select a story to preview its details here.

Module: Patient Registration
User Story ID: HC-089
User Story Title: Automated Unique Hospital Identification (UHID) Generation

User Story Description:
As a front-desk admission officer, I want the system to automatically generate a unique 10-digit UHID upon submitting new patient demographic details with barcode generation, so that every patient is distinctly trackable across all hospital departments (OPD, IPD, Diagnostics, Billing).

▶ Generate Without RAG (2 Cases)

▶ 🔧 Run Both Modes (Auto-Compare)

↔️ View Side-by-Side Comparison (Ready)

💾 Save Config

▶ Generated Test Cases Without RAG (Direct Baseline)

Mode: Without RAG Model: openai/gpt-oss-120b 0 DB Retrieval (Zero Context) ⚡ 2,556 / 8000 Tok (32% Budget Safe)

Export CSV

ID	Title	Module	Priority	Type	Preconditions	Test Steps	Expected Results
TC_6355	Positive: Successful UHID Generation with Valid Demographics	Patient Registration	P1	Functional	User is logged into the Hospital Management System with Admission Officer role; the Patient Registration screen is open; no existing patient record matches the provided national ID or phone number.	987654321V). 2. Fill mandatory fields: First Name, Last Name, Date of Birth, Phone Number (e.g., 555-123-4567), Emergency Contact Name/Phone, Blood Group (e.g., O+). 3. Click the "Submit" button. 4. Observe the system processing indicator. 5. Capture the generated UHID displayed on the confirmation screen. 6. Click the "Print Barcode" icon and verify that a barcode image appears. 7. Navigate to the patient profile and confirm that the UHID, barcode, and all entered data are persisted.	The system generates a UHID that matches the pattern <code>*[A-Z]{2}[d](8)\$</code> (two-letter prefix followed by eight digits) and the checksum digit is valid; no duplicate-national-ID warning appears; a 1-D barcode representing the UHID is rendered and is printable; patient record stores UHID, emergency contact, and blood group correctly; operation completes within 3 seconds.